

# **CQNTRM**

Datasheet revision 1.0 www.chipquik.com

# Flux Remover (Non-Flammable) in 10ml (0.34oz) Pen w/tip

## Product Highlights

A Flux Remover that is Non-Flammable Can be used with Leaded and Lead-Free applications RoHS II and REACH compliant



Packaging: 10ml Pen

Shelf Life: Refrigerated >24 months, Unrefrigerated >24 months



### Storage and Handling

Store refrigerated or at room temperature 3-25°C (37-77°F). The flux remover can be allowed to freeze. Freezing will not degrade this product. Allow 4 hours for flux remover to reach an operating temperature of 20-25°C (68-77°F) before use.

#### Transportation

This product has no shipping restrictions. Shipping below 0°C (32°F) or above 25°C (77°F) for normal transit times by ground or air will not impact this product's stated shelf life.

Conforms to the following Industry Standards:

RoHS 2 Directive 2011/65/EU:

Yes